



GGM

GGM HF-2900 Humidoflo® HFT System represents a new state of the art technology for delivery of high standard high flow therapy for pediatric and adult patients in the hospital, home, acute and subacute clinical settings.

Consists of a Servo Controlled Flow generator with an integrated humidification and oxygen enrichment systems.

Delivers up to 60 lpm of inspiratory flow in adult patients. This flow rate helps meet patients breathing requirements as well as help maintain their body temperature whilst delivering up to 99% relative humidity.

Small, Lightweight and intuitive, the **GGM HF-2900 Humidoflo® HFT System** is the ideal choice for the discerning buyer.



www.greatgroup.com.tw

HUMIDOFLO® HFT SYSTEM HF-2900



CE₂₄₆₀